

FIBP Antibody, HRP conjugated



PACO26774

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in mitogenic function of FGF1.

Gene ID:

FIBP

Uniprot

O43427

Synonyms:

acid, c fibroblast growth factor intracellular-binding protein (aFGF intracellular-binding protein) (FGF-1 intracellular-binding protein), FIBP

Immunogen:

Recombinant Human acid, c fibroblast growth factor intracellular-binding protein (2-357AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A