

# FIBP Antibody, Biotin conjugated



PACO26776

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in mitogenic function of FGF1.

**Gene ID:**

FIBP

**Uniprot**

O43427

**Synonyms:**

acid, c fibroblast growth factor intracellular-binding protein (aFGF intracellular-binding protein) (FGF-1 intracellular-binding protein), FIBP

**Immunogen:**

Recombinant Human acid, c fibroblast growth factor intracellular-binding protein (2-357AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A