

# FSCN2 Antibody, Biotin conjugated



PACO26840

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Acts as an actin bundling protein. May play a pivotal role in photoreceptor cell-specific events, such as disk morphogenesis.

**Gene ID:**

FSCN2

**Uniprot**

O14926

**Synonyms:**

Fascin-2 (Retinal fascin), FSCN2

**Immunogen:**

Recombinant Human Fascin-2 protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A