

# FXYP1 Antibody, Biotin conjugated



PACO26856

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May have a functional role in muscle contraction. Induces a hyperpolarization-activated chloride current when exogenously expressed.

**Gene ID:**

FXYP1

**Uniprot**

O00168

**Synonyms:**

Phospholemman (FXYP domain-containing ion transport regulator 1) (Sodium/potassium-transporting ATPase subunit FXYP1), FXYP1, PLM

**Immunogen:**

Recombinant Human Phospholemman protein (57-92AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A