

FXYP7 Antibody, FITC conjugated



PACO26859

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

integral component of plasma membrane, plasma membrane, sodium channel regulator activity, ion transmembrane transport, regulation of sodium ion transmembrane transporter activity, transmembrane transport.

Gene ID:

FXYP7

Uniprot

P58549

Synonyms:

FXYP domain-containing ion transport regulator 7, FXYP7

Immunogen:

Recombinant Human FXYP domain-containing ion transport regulator 7 protein (47-80AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A