

G0S2 Antibody, Biotin conjugated



PACO26868

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes apoptosis by binding to BCL2, hence preventing the formation of protective BCL2-BAX heterodimers.

Gene ID:

G0S2

Uniprot

P27469

Synonyms:

G0/G1 switch protein 2 (G0/G1 switch regulatory protein 2) (Putative lymphocyte G0/G1 switch gene), G0S2

Immunogen:

Recombinant Human G0/G1 switch protein 2 protein (1-103AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A