

# GAGE5 Antibody



PACO26901

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

GAGE5

**Uniprot**

Q13069

**Synonyms:**

G antigen 5 (GAGE-5) (Cancer/testis antigen 4.5) (CT4.5), GAGE5

**Immunogen:**


Recombinant Human G antigen 5 protein (1-117AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO26901 at dilution of 1:100.