

GALK2 Antibody, FITC conjugated



PACO26907

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts on GalNAc. Also acts as a galactokinase when galactose is present at high concentrations. May be involved in a salvage pathway for the reutilization of free GalNAc derived from the degradation of complex carbohydrates.

Gene ID:

GALK2

Uniprot

Q01415

Synonyms:

N-acetylgalactosamine kinase (EC 2.7.1.157) (GalNAc kinase) (Galactokinase 2), GALK2, GK2

Immunogen:

Recombinant Human N-acetylgalactosamine kinase protein (161-458AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A