

# GALK2 Antibody, Biotin conjugated



PACO26908

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Acts on GalNAc. Also acts as a galactokinase when galactose is present at high concentrations. May be involved in a salvage pathway for the reutilization of free GalNAc derived from the degradation of complex carbohydrates.

**Gene ID:**

GALK2

**Uniprot**

Q01415

**Synonyms:**

N-acetylgalactosamine kinase (EC 2.7.1.157) (GalNAc kinase) (Galactokinase 2), GALK2, GK2

**Immunogen:**

Recombinant Human N-acetylgalactosamine kinase protein (161-458AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A