

PACO26917

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

Gamt

Uniprot

P10868

Synonyms:

Guanidinoacetate N-methyltransferase (EC 2.1.1.2), Gamt


Immunogen:

Recombinant Rat Guanidinoacetate N-methyltransferase protein (1-236AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: Gamt antibody at 2 μ g/ml + Rat liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 27 kDa

Observed band size: 27 kDa