

# Gamt Antibody, HRP conjugated



PACO26918

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Gamt

**Uniprot**

P10868

**Synonyms:**

Guanidinoacetate N-methyltransferase (EC 2.1.1.2), Gamt

**Immunogen:**

Recombinant Rat Guanidinoacetate N-methyltransferase protein (1-236AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A