

GDAP1L1 Antibody, Biotin conjugated



PACO26936

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasm, glutathione transferase activity, glutathione metabolic process.

Gene ID:

GDAP1L1

Uniprot

Q96MZ0

Synonyms:

Ganglioside-induced differentiation-associated protein 1-like 1 (GDAP1-L1), GDAP1L1

Immunogen:

Recombinant Human Ganglioside-induced differentiation-associated protein 1-like 1 protein (101-260AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A