

Gdf5 Antibody, HRP conjugated



PACO26938

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Could be involved in bone and cartilage formation. Chondrogenic signaling is mediated by the high-affinity receptor BMPR1B.

Gene ID:

Gdf5

Uniprot

P43027

Synonyms:

Growth/differentiation factor 5 (GDF-5) (Bone morphogenetic protein 14) (BMP-14), Gdf5, Bmp14 Bp Gdf-5

Immunogen:

Recombinant Mouse Growth/differentiation factor 5 protein (376-495AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A