

TSTA3 Antibody, Biotin conjugated



PACO26952

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

TSTA3

Uniprot

Q13630

Synonyms:

GDP-L-fucose synthase (EC 1.1.1.271) (GDP-4-keto-6-deoxy-D-mannose-3,5-epimerase-4-reductase) (Protein FX) (Red cell NADP(H)-binding protein) (Short-chain dehydrogenase/reductase family 4E member 1), TSTA3, SDR4E1

Immunogen:

Recombinant Human GDP-L-fucose synthase protein (1-314AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A