

# GHRH Antibody



PACO26973

---

## Product Information

**Size:**

50ug

**Reactivity:**

Sheep

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

GRF is released by the hypothalamus and acts on the adenohypophyse to stimulate the secretion of growth hormone.

**Gene ID:**

GHRH

**Uniprot**

P07217

**Synonyms:**

Somatoliberin (Growth hormone-releasing factor) (GRF) (Growth hormone-releasing hormone), GHRH

**Immunogen:**

Recombinant Sheep Somatoliberin protein (1-44AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A