

GHRH Antibody, FITC conjugated



PACO26975

Product Information

Size:

50ug

Reactivity:

Sheep

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

GRF is released by the hypothalamus and acts on the adenohypophyse to stimulate the secretion of growth hormone.

Gene ID:

GHRH

Uniprot

P07217

Synonyms:

Somatoliberin (Growth hormone-releasing factor) (GRF) (Growth hormone-releasing hormone), GHRH

Immunogen:

Recombinant Sheep Somatoliberin protein (1-44AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A