

GJA4 Antibody, Biotin conjugated



PACO26988

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

Gene ID:

GJA4

Uniprot

P35212

Synonyms:

Gap junction alpha-4 protein (Connexin-37) (Cx37), GJA4

Immunogen:

Recombinant Human Gap junction alpha-4 protein (231-333AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A