

# GJB3 Antibody, HRP conjugated



PACO26990

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

**Gene ID:**

GJB3

**Uniprot**

O75712

**Synonyms:**

Gap junction beta-3 protein (Connexin-31) (Cx31), GJB3, CX31

**Immunogen:**

Recombinant Human Gap junction beta-3 protein (211-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A