

GK2 Antibody, Biotin conjugated



PACO27008

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Key enzyme in the regulation of glycerol uptake and metabolism.

Gene ID:

GK2

Uniprot

Q14410

Synonyms:

Glycerol kinase 2 (GK 2) (Glycerokinase 2) (EC 2.7.1.30) (ATP: glycerol 3-phosphotransferase 2) (Glycerol kinase, testis specific 2), GK2, GKP2 GKTA

Immunogen:

Recombinant Human Glycerol kinase 2 protein (2-553AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A