

GLP1R Antibody, Biotin conjugated



PACO27028

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This is a receptor for glucagon-like peptide 1. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.

Gene ID:

GLP1R

Uniprot

P43220

Synonyms:

Glucagon-like peptide 1 receptor (GLP-1 receptor) (GLP-1-R) (GLP-1R), GLP1R

Immunogen:

Recombinant Human Glucagon-like peptide 1 receptor protein (24-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A