

GLRX Antibody, FITC conjugated



PACO27035

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Has a glutathione-disulfide oxidoreductase activity in the presence of NADPH and glutathione reductase. Reduces low molecular weight disulfides and proteins.

Gene ID:

GLRX

Uniprot

P35754

Synonyms:

Glutaredoxin-1 (Thioltransferase-1) (TTase-1), GLRX, GRX

Immunogen:

Recombinant Human Glutaredoxin-1 protein (1-106AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A