

## Product Information

**Size:**

50ug

**Reactivity:**

Bovine

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Accelerates the intermembrane transfer of various glycolipids. Catalyzes the transfer of various glycosphingolipids between membranes but does not catalyze the transfer of phospholipids. May be involved in the intracellular translocation of glucosylceramides.

**Gene ID:**

GLTP

**Uniprot**

P68265

**Synonyms:**

Glycolipid transfer protein (GLTP), GLTP

**Immunogen:**

Recombinant Bovine Glycolipid transfer protein (2-209AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A