

GNAI3 Antibody, HRP conjugated



PACO27066

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems. G(k) is the stimulatory G protein of receptor-regulated K⁺ channels. The active GTP-bound form prevents the association of RGS14 with centrosomes and is required for the translocation of RGS14 from the cytoplasm to the plasma membrane. May play a role in cell division.

Gene ID:

GNAI3

Uniprot

P08754

Synonyms:

Guanine nucleotide-binding protein G(k) subunit alpha (G(i) alpha-3), GNAI3

Immunogen:

Recombinant Human Guanine nucleotide-binding protein G(k) subunit α protein (1-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A