

Gnao1 Antibody, Biotin conjugated



PACO27072

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Guanine nucleotide-binding proteins (G proteins) are involved as modulators or transducers in various transmembrane signaling systems. Stimulated by RGS14. The G(o) protein function is not clear.

Gene ID:

Gnao1

Uniprot

P18872

Synonyms:

Guanine nucleotide-binding protein G(o) subunit alpha, Gnao1, Gna0 Gnao

Immunogen:

Recombinant Mouse Guanine nucleotide-binding protein G(o) subunit α protein (2-354AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A