

GOSR2 Antibody



PACO27117

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200

Protein Background:

Involved in transport of proteins from the cis/medial-Golgi to the trans-Golgi network.

Gene ID:

GOSR2

Uniprot

O14653

Synonyms:

Golgi SNAP receptor complex member 2 (27 kDa Golgi SNARE protein) (Membrin),
GOSR2, GS27


Immunogen:

Recombinant Human Golgi SNAP receptor complex member 2 protein (1-190AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot

All lanes: GOSR2 antibody at 4 μ g/ml + Jurkat whole cell lysate


Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 25, 23 kDa

Observed band size: 25 kDa

.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO27117 at dilution of 1:100.