

# GPD2 Antibody, Biotin conjugated



PACO27148

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

mitochondrial inner membrane, calcium ion binding, glycerol-3-phosphate dehydrogenase activity, cellular lipid metabolic process, small molecule metabolic process.

**Gene ID:**

GPD2

**Uniprot**

P43304

**Synonyms:**

Glycerol-3-phosphate dehydrogenase, mitochondrial (GPD-M) (GPDH-M) (EC 1.1.5.3) (mtGPD), GPD2

**Immunogen:**

Recombinant Human Glycerol-3-phosphate dehydrogenase, mitochondrial protein (43-727AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A