

SNRPA Antibody, Biotin conjugated



PACO27160

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds stem loop II of U1 snRNA. It is the first snRNP to interact with pre-mRNA. This interaction is required for the subsequent binding of U2 snRNP and the U4/U6/U5 tri-snRNP. In a snRNP-free form (SF-A) may be involved in coupled pre-mRNA splicing and polyadenylation process.

Gene ID:

SNRPA

Uniprot

P09012

Synonyms:

U1 small nuclear ribonucleoprotein A (U1 snRNP A) (U1-A) (U1A), SNRPA

Immunogen:

Recombinant Human U1 small nuclear ribonucleoprotein A protein (8-278AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A