

Gstp1 Antibody, HRP conjugated



PACO27231

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Can metabolize 1-chloro-2,4-dinitrobenzene. Regulates negatively CDK5 activity via p25/p35 translocation to prevent neurodegeneration.

Gene ID:

Gstp1

Uniprot

P19157

Synonyms:

Glutathione S-transferase P 1 (Gst P1) (EC 2.5.1.18) (GST YF-YF) (GST class-pi) (GST-piB) (Preadipocyte growth factor), Gstp1, Gstpib

Immunogen:

Recombinant Mouse Glutathione S-transferase P 1 protein (2-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A