

# Gstp1 Antibody, Biotin conjugated



PACO27235

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Can metabolize 1-chloro-2,4-dinitrobenzene. Regulates negatively CDK5 activity via p25/p35 translocation to prevent neurodegeneration.

**Gene ID:**

Gstp1

**Uniprot**

P19157

**Synonyms:**

Glutathione S-transferase P 1 (Gst P1) (EC 2.5.1.18) (GST YF-YF) (GST class-pi) (GST-piB) (Preadipocyte growth factor), Gstp1, Gstpib

**Immunogen:**

Recombinant Mouse Glutathione S-transferase P 1 protein (2-210AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A