

Gstp1 Antibody, Biotin conjugated



PACO27236

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Regulates negatively CDK5 activity via p25/p35 translocation to prevent neurodegeneration.

Gene ID:

Gstp1

Uniprot

P04906

Synonyms:

Glutathione S-transferase P (EC 2.5.1.18) (Chain 7) (GST 7-7) (GST class-pi), Gstp1

Immunogen:

Recombinant Rat Glutathione S-transferase P protein (2-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A