

GTF2E2 Antibody, Biotin conjugated



PACO27252

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GTF2E2

Uniprot

P29084

Synonyms:

Transcription initiation factor IIE subunit beta (TFIIE-beta) (General transcription factor IIE subunit 2), GTF2E2, TF2E2

Immunogen:

Recombinant Human Transcription initiation factor IIE subunit β protein (1-291AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A