

HDGF Antibody, Biotin conjugated



PACO27272

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Heparin-binding protein, with mitogenic activity for fibroblasts. Acts as a transcriptional repressor.

Gene ID:

HDGF

Uniprot

P51858

Synonyms:

Hepatoma-derived growth factor (HDGF) (High mobility group protein 1-like 2) (HMG-1L2), HDGF, HMG1L2

Immunogen:

Recombinant Human Hepatoma-derived growth factor protein (1-239AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A