

PACO27277

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This protein is a minor sialoglycoprotein in human erythrocyte membranes.

Gene ID:

GYPE

Uniprot

P15421

Synonyms:

Glycophorin-E, GYPE, GPE

Immunogen:

Recombinant Human Glycophorin-E protein (20-52AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A