

# HBD Antibody, HRP conjugated



PACO27322

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in oxygen transport from the lung to the various peripheral tissues.

**Gene ID:**

HBD

**Uniprot**

P02042

**Synonyms:**

Hemoglobin subunit delta (Delta-globin) (Hemoglobin delta chain), HBD

**Immunogen:**

Recombinant Human Hemoglobin subunit delta protein (2-147AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A