

# HBQ1 Antibody



PACO27333

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

hemoglobin complex, oxygen transport.

**Gene ID:**

HBQ1

**Uniprot**

P09105

**Synonyms:**

Hemoglobin subunit theta-1 (Hemoglobin theta-1 chain) (Theta-1-globin), HBQ1

**Immunogen:**

Recombinant Human Hemoglobin subunit theta-1 protein (1-142AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A