

HEMK1 Antibody, HRP conjugated



PACO27346

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

N5-glutamine methyltransferase responsible for the methylation of the GGQ triplet of the mitochondrial translation release factor MTRF1L.

Gene ID:

HEMK1

Uniprot

Q9Y5R4

Synonyms:

HemK methyltransferase family member 1 (EC 2.1.1. -) (M. HsaHemKP), HEMK1, HEMK

Immunogen:

Recombinant Human HemK methyltransferase family member 1 protein (1-338AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A