

HEXIM2 Antibody, HRP conjugated



PACO27354

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional regulator which functions as a general RNA polymerase II transcription inhibitor. In cooperation with 7SK snRNA sequesters P-TEFb in a large inactive 7SK snRNP complex preventing RNA polymerase II phosphorylation and subsequent transcriptional elongation.

Gene ID:

HEXIM2

Uniprot

Q96MH2

Synonyms:

Protein HEXIM2 (Hexamethylene bis-acetamide-inducible protein 2), HEXIM2

Immunogen:

Recombinant Human Protein HEXIM2 protein (1-286AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A