

# DAP Antibody, Biotin conjugated



PACO27360

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Involved in mediating interferon-gamma-induced cell death.

**Gene ID:**

DAP

**Uniprot**

P51397

**Synonyms:**

Death-associated protein 1 (DAP-1), DAP, DAP1

**Immunogen:**

Recombinant Human Death-associated protein 1 protein (2-102AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A