

HNRNPL Antibody, HRP conjugated



PACO27414

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Splicing factor binding to exonic or intronic sites and acting as either an activator or repressor of exon inclusion. Exhibits a binding preference for CA-rich elements. Component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complexes and associated with most nascent transcripts. Associates, together with APEX1, to the negative calcium responsive element (nCaRE) B2 of the APEX2 promoter.

Gene ID:

HNRNPL

Uniprot

P14866

Synonyms:

Heterogeneous nuclear ribonucleoprotein L (hnRNP L), HNRNPL, HNRPL

Immunogen:

Recombinant Human Heterogeneous nuclear ribonucleoprotein L protein (89-335AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A