

# Hpx Antibody, Biotin conjugated



PACO27444

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Binds heme and transports it to the liver for breakdown and iron recovery, after which the free hemopexin returns to the circulation.

**Gene ID:**

Hpx

**Uniprot**

P20059

**Synonyms:**

Hemopexin, Hpx

**Immunogen:**

Recombinant Rat Hemopexin protein (24-460AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A