

IFNA21 Antibody, HRP conjugated



PACO27534

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	IFNA21
Source:	Uniprot
Rabbit	P01568
Isotype:	Synonyms:
IgG	Interferon alpha-21 (IFN-alpha-21) (Interferon alpha-F) (LeIF F), IFNA21
Applications:	Immunogen:
ELISA	Recombinant Human Interferon alpha-21 protein (24-189AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A