

TTC1 Antibody, FITC conjugated



PACO27543

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Interacts with the GAP domain of NF1.

Gene ID:

TTC1

Uniprot

Q99614

Synonyms:

Tetratricopeptide repeat protein 1 (TPR repeat protein 1), TTC1, TPR1

Immunogen:

Recombinant Human Tetratricopeptide repeat protein 1 protein (3-292AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A