

IGSF6 Antibody, HRP conjugated



PACO27574

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

integral component of plasma membrane, transmembrane signaling receptor activity, cell surface receptor signaling pathway, immune response.

Gene ID:

IGSF6

Uniprot

O95976

Synonyms:

Immunoglobulin superfamily member 6 (IgSF6) (Protein DORA), IGSF6, DORA

Immunogen:

Recombinant Human Immunoglobulin superfamily member 6 protein (55-153AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A