

IL2RG Antibody, Biotin conjugated



PACO27616

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Common subunit for the receptors for a variety of interleukins.

Gene ID:

IL2RG

Uniprot

P31785

Synonyms:

Cytokine receptor common subunit gamma (Interleukin-2 receptor subunit gamma) (IL-2 receptor subunit gamma) (IL-2R subunit gamma) (IL-2RG) (gammaC) (p64) (CD antigen CD132), IL2RG

Immunogen:

Recombinant Human Cytokine receptor common subunit γ protein (75-230AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A