

PACO27669

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Transcription factor that regulates gene expression in intestine. May participate in vitamin A metabolism most likely by regulating BCO1 expression in the intestine.

**Gene ID:**

ISX

**Uniprot**

Q2M1V0

**Synonyms:**

Intestine-specific homeobox (RAX-like homeobox), ISX, RAXLX

**Immunogen:**


Recombinant Human Intestine-specific homeobox protein (1-245AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO27669 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human glioma using PACO27669 at dilution of 1:100.