

# JAK2 Antibody, FITC conjugated



PACO27715

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

JAK2

**Uniprot**

O60674

**Synonyms:**

Tyrosine-protein kinase JAK2 (EC 2.7.10.2) (Janus kinase 2) (JAK-2), JAK2

**Immunogen:**

Recombinant Human Tyrosine-protein kinase JAK2 protein (752-1132AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A