

JAM2 Antibody, Biotin conjugated



PACO27724

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in the processes of lymphocyte homing to secondary lymphoid organs.

Gene ID:

JAM2

Uniprot

P57087

Synonyms:

Junctional adhesion molecule B (JAM-B) (Junctional adhesion molecule 2) (JAM-2) (Vascular endothelial junction-associated molecule) (VE-JAM) (CD antigen CD322), JAM2, C21orf43 VEJAM

Immunogen:

Recombinant Human Junctional adhesion molecule B protein (29-238AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A