

# IMPDH2 Antibody, FITC conjugated



PACO27727

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:****Gene ID:**

IMPDH2

**Uniprot**

P12268

**Synonyms:**

Inosine-5'-monophosphate dehydrogenase 2 (IMP dehydrogenase 2) (IMPD 2) (IMPDH 2) (EC 1.1.1.205) (IMPDH-II), IMPDH2, IMPD2

**Immunogen:**

Recombinant Human Inosine-5'-monophosphate dehydrogenase 2 protein (5-514AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A