

KCNJ1 Antibody, Biotin conjugated



PACO27752

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

KCNJ1

Uniprot

P48048

Synonyms:

ATP-sensitive inward rectifier potassium channel 1 (ATP-regulated potassium channel ROM-K) (Inward rectifier K(+) channel Kir1.1) (Potassium channel, inwardly rectifying subfamily J member 1), KCNJ1, ROMK1

Immunogen:

Recombinant Human ATP-sensitive inward rectifier potassium channel 1 protein (178-391AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A