

KIAA1324L Antibody



PACO27789

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

KIAA1324L

Uniprot

A8MWY0

Synonyms:

UPF0577 protein KIAA1324-like (Estrogen-induced gene 121-like protein) (hEIG121L), KIAA1324L, EIG121L


Immunogen:

Recombinant Human UPF0577 protein KIAA1324-like protein (227-526AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO27789 at dilution of 1:100.